Head-on PMT

Photon Counting Head H9319 Series



The H9319 series photon counting heads are designed to perform photon counting by just connecting to a PC. The H9319 series includes a 25-mm (1") diameter head-on photomultiplier tube, a photon counting circuit, a high-voltage power supply circuit, counter and a microprocessor. Data transfer, measurement time and other necessary adjustments can be controlled by commands from the PC through the RS-232C interface.

Since the H9319 series performs linearity correction by the internal microprocessor, it provides excellent count linearity within a range of ± 1 % at 20×10^6 s⁻¹.

Product Variations

Type No.	Spectral Response	Sample Program	Prescaler
H9319-01	300 nm to 650 nm	yes	1/4
H9319-11		no	1/4
H9319-02	300 nm to 850 nm	yes	1/4
H9319-12		no	1/4

This product can't be used at vacuum environment or reduced pressure environment.

Specifications

(at +25 °C)

Parameter			H9319-01	H9319-11	H9319-02	H9319-12	Unit	
Input Voltage		+4.75 to +5.25				V		
Max. Input Voltage			+6			V		
Max. Input Current			60			mA		
Effective Area			φ 22			mm		
Peak Sensitivity Wavelength		420			nm			
Count Sensitivity	Тур.	300 nm	2.3 >	< 10 ⁵	2.1 >	< 10 ⁵		
		400 nm	4.1×10^5 $2.5 \times$		< 10 ⁵			
		500 nm	3.4×10^5 2.0×10^5		s-¹⋅pW-¹			
		600 nm	5.7 × 10 ⁴		1.3 × 10 ⁵			
		700 nm	_		7.8 >	< 10 ⁴		
Count Linearity *1		20 × 10 ⁶			S ⁻¹			
Dark Count *2		Тур.	150		10 000		- s ⁻¹	
		Max.	300		15 000			
PMT Operating Voltage Range		+300 to +1200			V			
Integration Time		10 to 1000			ms			
Settling Time		1 *3			S			
		5 *4			S			
Input Signal (External Trigger Input) *5		TTL level signal			_			
Output Signal (User Line Output) *6		TTL level signal			_			
Interface			RS-232C, 9600 baud, Parity none, 8 data bit, 1 stop bit			_		
Operating Ambient Temperature *7			+5 to +50			°C		
Storage Temperature *7			-20 to +50			°C		
Weight *8		280			g			

^{*1:} Random pulse, within ±1 % count loss (by count linearity compensation)

^{*2:} After 30 minutes storage in darkness

^{*3:} The time required for the output to reach a stable level following a change in the control voltage from 500 V to 1000 V in darkness

^{*4:} The time required for the output to reach a stable level following a change in the control voltage from 1000 V to 500 V in darkness

^{*5:} Selectable ECR mode or LCR mode (refer to command list)


^{*6:} Controllable by RS-232C command

^{*7:} No condensation


^{*8:} Main body: Approx. 180 g

Photon Counting Head with CPU+Interface

Characteristic (Count sensitivity)



Block Diagram




Command List

Action	Command *9	Explanation	Argument	Response *10
Set the Integration Time	P#C _R	Set the number of 10 msec intervals to sum. Same as integration time for 1 reading.	1 to 100	VA, BC, BA
Set the Sequence		Set sequence of readings, where each reading uses the integration time set with the P command.	1 to 255	VA, BC, BA
Set the PMT Input	V##CR	Change the high voltage applied to the tube.	300 to 1200	VA, BC, BA
Voltage	DCR	Re-set the default high voltage to the tube	_	VA, BC
Set the Output of the User Line	O#CR	Set the output of the user digital output line. 0: Low level output 1: High level output	0 or 1	VA, BC, BA
	SCR	Start the reading sequence	_	4 byte/reading
Chart the Deading	CCR Start a continuous reading process. Will continue indefinitely until a STOP character (CR) is sent.		_	4 byte/reading
Start the Reading Sequence	ECR	Start a reading sequence for each positive-edge TTL transition applied to the purple stripe user line.	_	4 byte/reading
	LCR	Start a reading sequence for each positive-level TTL transition applied to the purple stripe user line.	_	4 byte/reading

^{*9:} The Response acknowledgment is returned having two bytes.

Dimensional Outlines (Unit: mm)



65

^{*10:} VA: valid command, BC: bad command, BA: bad argument